

Rabbit anti-Estriol

| | | | | | | | | | | | | | |
|--|---|---------|------|-----------|------|---------------|------|--|------|---------|-------|----------------|-------|
| CATALOG NUMBER: | P124 | | | | | | | | | | | | |
| HOST ANIMAL: | Rabbit | | | | | | | | | | | | |
| SOURCE: | Serum from rabbits immunized with Estriol-BSA | | | | | | | | | | | | |
| FORM: | Neat antiserum. Contains 0.1% NaN ₃ | | | | | | | | | | | | |
| CROSS REACTIVITY: | <table><tr><td>Estriol</td><td>100%</td></tr><tr><td>Estradiol</td><td>0.1%</td></tr><tr><td>16-Epiestriol</td><td>0.2%</td></tr><tr><td>Estriol-3-H₂SO₄</td><td>0.3%</td></tr><tr><td>Estrone</td><td>0.10%</td></tr><tr><td>Other steroids</td><td><0.1%</td></tr></table> | Estriol | 100% | Estradiol | 0.1% | 16-Epiestriol | 0.2% | Estriol-3-H ₂ SO ₄ | 0.3% | Estrone | 0.10% | Other steroids | <0.1% |
| Estriol | 100% | | | | | | | | | | | | |
| Estradiol | 0.1% | | | | | | | | | | | | |
| 16-Epiestriol | 0.2% | | | | | | | | | | | | |
| Estriol-3-H ₂ SO ₄ | 0.3% | | | | | | | | | | | | |
| Estrone | 0.10% | | | | | | | | | | | | |
| Other steroids | <0.1% | | | | | | | | | | | | |
| PRESERVATIVE: | 0.1% sodium azide | | | | | | | | | | | | |
| STORAGE: | Store antibody at 2-8°C short term and -20°C long term | | | | | | | | | | | | |
| INTENDED USE: | For Research Use Only. Not intended for diagnostic use. Suitable for use in immunoassay development or other applications. | | | | | | | | | | | | |